



Affordable | Profitable | Reliable | Effective | User Friendly

Sapphire

LAMP PUMPED SERIES

Sapphire
YAG



THE SAPPHIRE Nd:YAG 1064

The ideal solution for the busy clinic, which can be safely and successfully used for the removal of hair on skin types I-VI

Laser light at 1064nm penetrates deeper into the skin enabling the targeting of a broader range of skin and hair types, making long pulsed Nd: Yag lasers the best choice for hair removal on skin types I-VI.

Unwanted hair can be a thing of the past using the Sapphire long pulsed Nd: Yag Laser.

Permanent hair reduction can be successfully achieved for both men and women.

All areas of the body can be safely treated, using the ergonomic and easy to use hand piece.

Made in the U.K.

SPECIFICATION

| | |
|-----------------|-------------------------|
| Wavelength | 1064nm |
| Mode | Long pulsed |
| Energy | up to 45J per pulse |
| Pulse duration | 10ms |
| Repetition rate | up to 1.7Hz |
| Beam delivery | Fibre optic |
| Spot size | 4mm & 10mm |
| Aiming Beam | 635nm Diode |
| Fluence | upto 300 J/cm sq |
| Size | H 88, W 35, D 78cm |
| Weight | 65Kg |
| Electrical | 240V 113A / 50HZ |
| System Cooling | Internal self contained |



Sapphire
ALEX



THE SAPPHIRE ALEXANDRITE

The Sapphire Alexandrite Laser remains one of the best laser tools for the removal of hair for skin types I-III

With high power per pulse and variable spot sizes treatments are quick and effective.

Unwanted hair can be a thing of the past using the Sapphire long pulsed Alexandrite laser.

Permanent hair reduction can be successfully achieved for both men and women.

All areas of the body can be safely treated, using the ergonomic and easy to use hand piece.

Made in the U.K.

SPECIFICATION

| | |
|-----------------|-----------------------------|
| Wavelength | 755nm |
| Mode | Long pulsed |
| Energy | up to 50J per pulse |
| Pulse duration | 3ms |
| Repetition rate | 1Hz |
| Beam delivery | Fibre optic |
| Spot size | 10mm, 15mm & 18mm |
| Aiming Beam | 532nm Diode |
| Fluence | upto 30 J/cm sq (15mm Spot) |
| Size | H 88, W 40, D 80cm |
| Weight | 65 Kg |
| Electrical | 240V / 13A / 50HZ |
| System Cooling | Internal self contained |



Sapphire RUBY



THE SAPPHIRE RUBY 694

The Sapphire RUBY laser remain one of the best tools for the removal of tattoos and pigmented lesions.

With high powers and adjustable spot sizes, treatments are quick and effective.

The Sapphire Q- Switched Ruby Laser produces pulses of energy that are very short, but high in power.

The light from the laser is absorbed by the pigment in the tattoo but not the surrounding skin.

A combination of light absorption and acoustic shock help to break down the ink into smaller particles that the body can dispose of naturally.

Made in the U.K.

SPECIFICATION

| | |
|-----------------|-------------------------|
| Wavelength | 694nm |
| Mode | Q-Switched |
| Energy | up to 1.25J per pulse |
| Pulse duration | 35ns |
| Repetition rate | 1Hz |
| Beam delivery | Articulated arm |
| Spot size | 4mm, 5mm & 6mm |
| Aiming Beam | 633nm HeNe |
| Fluence | 1,8 to 9.95 J/cm sq |
| Size | H 102, W 40, D 93cm |
| Weight | 115 Kg |
| Electrical | 240V / 13A / 50HZ |
| System Cooling | Internal self contained |



A World Class Service from a World Class Company



Sapphire
ALEX

Remains one of the best
Laser Tools for removal
of hair for skin types I-III



Sapphire
RUBY

Remains one of the best tools
for the removal of tattoos and
pigmented lesions



Sapphire
YAG

The ideal solution for the
busy clinic which can
safely & successfully be
used for removal of hair
on skin types I-VI

LASER S.O.S. AESTHETICS

Laser SOS Aesthetics Ltd

Unit 3 Burrell Road, St. Ives, Cambs PE27 3LE, United Kingdom

Tel: +44 (0) 1480 409962

Fax: +44 (0) 1480 498860

Email: sales@lasersosaesthetics.com

Web: www.lasersosaesthetics.com